



Large Darkroom (B-L)



A darkroom shields light.
Then, it makes possible to stabilize of optical measurement.
The treatment of preventing dust and static electricity are given to our darkroom.
You can darken your work area locally.

For the wide field

Working with microscope, a safety measure against laser machine, Chemistry and the field of the material research ...etc



Make up easily by frame pipes!!

Time: 30min~1h
(Depends on the size)

Tool set

(Wrench, Screwdriver)

**Instruction manual
with pictures**

Custom made service

(Within 1.6~4.2meter)

List

We have three kinds of darkroom. Please select from following list.

**Hook-and-Loop
fastener type**



Curtain type



Cross type



This curtain is crossed in the middle.

Hook-and-Loop Fastener type	Curtain type	Cross type	Size (mm)
B-L1	B-L1-CU	B-L1-X	(W)1600×(D)1600×(H)2000
B-L2	B-L2-CU	B-L2-X	(W)2100×(D)1600×(H)2000
B-L3	B-L3-CU	B-L3-X	(W)2100×(D)2100×(H)2000
B-L4	B-L4-CU	B-L4-X	(W)3100×(D)1600×(H)2000
B-L5	B-L5-CU	B-L5-X	(W)3100×(D)2100×(H)2000
B-L6	B-L6-CU	B-L6-X	(W)4200×(D)1600×(H)2000
B-L7	B-L7-CU	B-L7-X	(W)4200×(D)2100×(H)2000
<u>Available at clean room</u>			We can provide special size on request. Width: 1.6~4.2meter / Depth 1.6-2.1m Height: 2m

Specification

<Frame>diameter 30mm,plastic coating iron pipe <weight>6-10kg(the total of the frame and the dark curtain)
 <The quality of dark curtain> The shield light ability is the first rank (100-99.99% shielding . Refer to a note 1 for the details.) <Front side> Polyester 100% Non-flammable material. <Rear side>Polyurethane black. Laminating processing. The treatment of preventing dust and static electricity.
 <Shielding light ability>At Visible ray, Ultraviolet ray and Infrared ray : more than 99.99%.(Refer to a note 1)
 *The skirt of dark curtain is made a little long to prevent optical omission.
 *Frame pipe and dark curtain are taken apart to ship them. The tools and the manual are packed with them.

*The specification of our product may be changed for the reason such as improvement and the progress of technology without notice.

2016/5/30

SCIENTEX Inc.

Sales department (Business hour 9:00~18:00)

1085-1 nishimisono hamakita-ku Hamamatsu-city Shizuoka-pref 434-0033 Japan

URL <http://www.scientex.co.jp> E-mail info@scientex.co.jp

Overseas agency